





## PRODUCT DATASHEET- STERILE LAMINATED SURGICAL GOWN

<b>Reference No</b>	40.00X.01
<b>Product Name</b>	Sterile Laminated Surgical Gown
<b>Trademark</b>	Medcare
<b>Manufacturer</b>	Medcare Sağlık Ürünleri San. Ve Tic. A.Ş.

Front side of laminated gown	Back side of laminated gown
 A blue, long-sleeved, full-length laminated surgical gown shown from the front. It has a V-neckline and a simple, smooth surface.	 A blue, long-sleeved, full-length laminated surgical gown shown from the back. It features a V-neckline and a large, central, rectangular opening at the waist, secured by a wide, flat fabric tie.

**PRODUCT DATA SHEET**  
**STERILE LAMINATED SURGICAL GOWN**  
**40.00X.01**

**1) PRODUCT COMPONENTS**

	
<p align="center"><b>ITEM:1</b>  <b>Sterile Laminated Surgical Gown</b>  <b>Size: S-M-L-XL-XXL-XXXL</b></p>	<p align="center"><b>ITEM: 2</b>  <b>Hand Towel</b>  <b>30x40 cm</b></p>

**2) DETAILS OF PRODUCT COMPONENTS**

QTY	PRODUCT	FABRIC TYPE	DIMENSION	COLOUR
1	Sterile Laminated Surgical Gown	Two Layer Fabric Biflex 56 gr/m <sup>2</sup> (+/-3) (30 gr PP Nonwoven+2 gr Hotmelt+24 gr PE Film)	S-M-L-XL-XXL-XXXL	Blue
1	Hand Towel	Cellulose	30x40 cm	White

**LAMINATED SURGICAL GOWN MEASUREMENT (cm)**

SIZE	CODE	A	B	C	D	E
<b>S</b>	40.001.01	110	100	28	55	140
<b>M</b>	40.002.01	115	105	29	58	145
<b>L</b>	40.003.01	127	117	29	60	150
<b>XL</b>	40.004.01	135	125	32	62	150
<b>XXL</b>	40.005.01	140	130	32	64	155
<b>XXXL</b>	40.006.01	145	135	34	66	160

Tolerance +/-%3 (cm)




**PRODUCT DATA SHEET**  
**STERILE LAMINATED SURGICAL GOWN**  
**40.00X.01**

**3) GENERAL INFORMATION**

<b>Product Description</b>	STERILE STANDARD SURGICAL GOWN
<b>Packing Quantity</b>	40 pcs
<b>Sterilization</b>	EO (Ethylene Oxide)
<b>Certificates</b>	EC Certificate
	ISO 13485:2016
<b>Standards</b>	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
<b>Instructions For Usage</b>	This product can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
<b>Expiry Date</b>	3 years

**PRODUCT DATA SHEET**  
**STERILE LAMINATED SURGICAL GOWN**  
**40.00X.01**

**4) PACKAGING**

<b>All production</b>	
Labelling of product	Clean Room Conditions, Medical Folding, Sterilization Bag, Peelable label,
Each unit	
<ul style="list-style-type: none"> <li>- Article number,</li> <li>- Lot number,</li> <li>- Manufacturer ,</li> <li>- Production date,</li> <li>- Expiry date,</li> <li>- EO Sterilization,</li> <li>- Single use,</li> <li>- CE,</li> <li>- Other customer requests,</li> </ul>	
<b>Packaging</b>	
<u>Pallet system</u>	<u>Carton</u>
Height (cm) -A	40
Width (cm) - B	60
Length (cm) - C	40
	
<b>Storage Instructions</b>	<p><b>HANDLING AND STORAGE:</b>            Products are sterile and must be stored in its original closed carton box, in a dry place. The products must be used as soon as the package is opened in the operation room. Do not expose to humidity and direct heat or sunlight.</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>% 10 Humidity Limitation</p> </div> </div>