



medcare

PRODUCT DATASHEET-STERILE REINFORCED SURGICAL GOWN

Reference No	20.00X.01
Product Name	Sterile Reinforced Surgical Gown
Trademark	Medcare
Manufacturer	Medcare Sağlık Ürünleri San. Ve Tic. A.Ş.

Front side of reinforced surgical gown	Back side of reinforced surgical gown
	

PRODUCT DATA SHEET
STERILE REINFORCED SURGICAL GOWN
 20.00X.01

1) PRODUCT COMPONENTS	
	
ITEM:1 Sterile Reinforced Surgical Gown Size: S-M-L-XL-XXL-XXXL	ITEM: 2 Hand Towel 30x40 cm

2) DETAILS OF PRODUCT COMPONENTS				
QTY	PRODUCT	PRODUCTION TYPE	DIMENSION	COLOUR
1	Sterile Reinforced Surgical Gown	Sms 43gr/m2 (+/-3) Reinforcement fabric 26gr/m2 (+/-3) (Spunbond PP/PE film)	Size: S-M-L-XL-XXL-XXXL	Blue
1	Hand Towel	Cellulose	30x40 cm	White

REINFORCED SURGICAL GOWN MEASUREMENT (cm)							
SIZE	A	B	C	D	E	F	G
S	110	100	28	55	140	40	60
M	115	105	29	58	145	45	65
L	127	117	29	60	150	50	70
XL	135	125	32	62	155	50	90
XXL	140	130	32	64	155	50	100
XXXL	145	135	33	66	160	55	110




PRODUCT DATA SHEET
STERILE REINFORCED SURGICAL GOWN
20.00X.01

3) GENERAL INFORMATION

Product Description	STERILE REINFORCED SURGICAL GOWN
Packing Quantity	S-M-L-XL 40 pcs XXL 35 pcs XXXL 30 pcs
Sterilization	EO (Ethylene Oxide)
Certificates	EC Certificate
	ISO 13485:2016
Standards	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
Instructions For Usage	This product can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
Expiry Date	3 years

PRODUCT DATA SHEET
STERILE REINFORCED SURGICAL GOWN
20.00X.01

4) PACKAGING

All production				
Labelling of product	Clean Room Conditions,	Medical Folding,	Sterilization Bag,	Peelable label,
Each unit <ul style="list-style-type: none"> - Article number, - Lot number, - Manufacturer , - Production date, - Expiry date, - EO Sterilization, - Single use, - CE, - Other customer requests, 				
Packaging				
<u>Pallet system</u>	<u>Carton</u>			
Height (cm) - A	40			
Width (cm) - B	60			
Length (cm) - C	40			
Storage Instructions	<p>HANDLING AND STORAGE: Products are sterile and must be stored in its original closed carton box, in a dry place. The products must be used as soon as the package is opened in the operation room. Do not expose to humidity and direct heat or sunlight.</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>Humidity Limitation</p> </div> </div>			