


**DELIVERY DRAPE PACK**  
**REF.NO:160.001.01**  
**STANDARD PRODUCT DATA SHEET**  
**KY-PD-030**



STANDARD PRODUCT DATA SHEET  
DELIVERY DRAPE PACK  
160.001.01


1) PRODUCT COMPONENTS



**ITEM 1**  
Fenestrated Drape  
130x150 cm




**ITEM 2**  
Leggings  
50x120 cm (unfolded)




**ITEM 3**  
Plain Drape  
100x100 cm


STANDARD PRODUCT DATA SHEET  
DELIVERY DRAPE PACK  
160.001.01




**ITEM 4**  
Baby Drape  
60x80 cm



**ITEM 5**  
Surgical Gown  
L



**ITEM 6**  
Hand Towel  
30x40 cm



**ITEM 7**  
Back Table Cover  
100x150 cm  
*(All of the product components are wrapped into back table cover)*

**STANDARD PRODUCT DATA SHEET  
DELIVERY DRAPE PACK  
160.001.01**

**2) DETAILS OF PRODUCT COMPONENTS**

ITEM	QTY	PRODUCT	FABRIC TYPE	DIMENSION	COLOUR
1	1	Fenestrated Drape	Biflex 56gr/m2 (+/-3)	130x150 cm	Blue
2	2	Leggings	Biflex 56gr/m2 (+/-3)	50x120 cm	Blue
3	1	Plain Drape	Biflex 56gr/m2 (+/-3)	100x100 cm	Blue
4	1	Baby Drape	Viscose	60x80 cm	White
5	2	Surgical Gown	Sms 43gr/m2 (+/-3)	L	Blue
6	2	Hand Towel	Cellulose	30x40 cm	White
7	1	Back Table Cover	Biflex 56gr/m2 (+/-3)	100x150 cm	Blue

**STANDARD SURGICAL GOWN MEASUREMENT (cm)**

SIZE	A	B	C	D	E
L	127	117	29	60	150

**STANDARD PRODUCT DATA SHEET  
DELIVERY DRAPE PACK  
160.001.01**

**3)GENERAL INFORMATION**

<b>Product Description</b>	<b>DELIVERY DRAPE PACK</b>
<b>Reference Code</b>	160.001.01
<b>Packing Quantity</b>	1 pc
<b>Quantity Per Box</b>	10 pcs
<b>Sterilization</b>	EO (Ethylene Oxide)
<b>Certificates</b>	EC Certificate
	ISO 13485:2016
<b>Standards</b>	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
<b>Instructions For Usage</b>	This set can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
<b>Expiry Date</b>	3 years

**STANDARD PRODUCT DATA SHEET**  
**DELIVERY DRAPE PACK**  
**160.001.01**

**5) PACKAGING**

<b>All production</b>	Clean Room Conditions, Medical Folding, Sterilization Bag, Peelable label,	
<b>Labelling of product</b>	Each unit <ul style="list-style-type: none"> <li>- Article number,</li> <li>- Lot number,</li> <li>- Manufacturer ,</li> <li>- Production date,</li> <li>- Expiry date,</li> <li>- EO Sterilization,</li> <li>- Single use,</li> <li>- CE,</li> <li>- Other customer requests,</li> </ul>	
<b>Packaging</b>	<u>Pallet system</u>	<u>Carton</u>
	Height (cm) -A	40
	Width (cm) - B	60
	Length (cm) - C	40
		
<b>Storage Instructios</b>	<p><b>HANDLING AND STORAGE:</b>  Products are sterile and must be stored in its original closed carton box, in a dry place.  The products must be used as soon as the package is opened in the operation room.  Do not expose to humidity and direct heat or sunlight.</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>Humidity Limitation</p> </div> </div>	