


**GENERAL SURGERY DRAPE PACK  
REFERENCE CODE: 110.001.01  
STANDARD PRODUCT DATA SHEET  
KY-PD-001**




STANDARD PRODUCT DATA SHEET  
GENERAL SURGERY DRAPE PACK  
110.001.01


1) PRODUCT COMPONENTS




**ITEM 1**  
Adhesive Drape  
(60 mm x 95 cm adhesive with fingerlift)  
150x240 cm



**ITEM 2**  
Adhesive Drape  
(60 mm x 95 cm adhesive with fingerlift)  
180x175 cm




**ITEM 3**  
Adhesive Drape  
(60 mm x 90 cm adhesive with fingerlift)  
75x90 cm




**ITEM 4**  
Mayo Stand Cover  
78x145 cm (unfolded)


**STANDARD PRODUCT DATA SHEET**  
**GENERAL SURGERY DRAPE PACK**  
**110.001.01**



**ITEM 5**  
Hand Towels  
30x40 cm



**ITEM 6**  
Op-Tape  
10x50 cm



**ITEM 7**  
Back Table Cover  
150x190 cm  
*(All of the product components are wrapped into back table cover)*

**STANDARD PRODUCT DATA SHEET**  
**GENERAL SURGERY DRAPE PACK**  
110.001.01

**2) DETAILS OF PRODUCT COMPONENTS**

ITEM	QTY	PRODUCT	FABRIC TYPE	DIMENSION	COLOUR
1	1	Adhesive Drape <i>(60mm x 95 cm adhesive with filgerlift)</i>	Biflex 56gr/m <sup>2</sup> (+/-3)	150x240 cm	Blue
2	1	Adhesive Drape <i>(60 mm x 95 cm adhesive with filgerlift)</i>	Biflex 56gr/m <sup>2</sup> (+/-3)	180x175 cm	Blue
3	2	Adhesive Drape <i>(60 mm x 90 cm adhesive with fingerlift)</i>	Biflex 56gr/m <sup>2</sup> (+/-3)	75x90 cm	Blue
4	1	Mayo Stand Cover	Pe+nonwoven	78x145 cm	Blue
5	4	Hand Towels	Cellulose	30x40 cm	White
6	1	OP-Tape	Sms+hotmelt+silicon paper	10x50 cm	Blue
7	1	Back Table Cover	Biflex 56gr/m <sup>2</sup> (+/-3)	150x190 cm	Blue
<i>All of the product components are wrapped into back table cover.</i>					




**STANDARD PRODUCT DATA SHEET**  
**GENERAL SURGERY DRAPE PACK**  
110.001.01

**3)GENERAL INFORMATION**

<b>Product Description</b>	<b>GENERAL SURGERY DRAPE PACK</b>
<b>Reference Code</b>	110.001.01
<b>Packing Quantity</b>	1 pc
<b>Quantity Per Box</b>	12 pcs
<b>Sterilization</b>	EO (Ethylene Oxide)
<b>Certificates</b>	EC Certificate
	ISO 13485:2016
<b>Standards</b>	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
<b>Instructions For Usage</b>	This set can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
<b>Expiry Date</b>	3 years

**STANDARD PRODUCT DATA SHEET**  
**GENERAL SURGERY DRAPE PACK**  
110.001.01

**5) PACKAGING**

<b>All production</b>	Clean Room Conditions, Medical Folding, Sterilization Bag, Peelable label,	
<b>Labelling of product</b>	Each unit <ul style="list-style-type: none"> <li>- Article number,</li> <li>- Lot number,</li> <li>- Manufacturer ,</li> <li>- Production date,</li> <li>- Expiry date,</li> <li>- EO Sterilization,</li> <li>- Single use,</li> <li>- CE,</li> <li>- Other customer requests,</li> </ul>	
<b>Packaging</b>	<u>Pallet system</u>	<u>Carton</u>
	Height (cm) - A	40
	Width (cm) - B	60
	Length (cm) - C	40
		
<b>Storage Instructios</b>	<p><b>HANDLING AND STORAGE:</b>  Products are sterile and must be stored in its original closed carton box, in a dry place.  The products must be used as soon as the package is opened in the operation room.  Do not expose to humidity and direct heat or sunlight.</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>Humidity Limitation</p> </div> </div>	