


**GENERAL SURGERY DRAPE PACK
REFERENCE CODE: 111.001.01
STANDARD PRODUCT DATA SHEET
KY-PD-002**




STANDARD PRODUCT DATA SHEET
GENERAL SURGERY DRAPE PACK
111.001.01


1) PRODUCT COMPONENTS



ITEM 1
Adhesive Drape
(60 mm x 95 cm adhesive with fingerlift)
150x240 cm



ITEM 2
Adhesive Drape
(60 mm x 95 cm adhesive with fingerlift)
150x180 cm





ITEM 3
Adhesive Drape
(60 mm x 95 cm adhesive with fingerlift)
75x90 cm




ITEM 4
Mayo Stand Cover
78x145 cm (unfolded)

STANDARD PRODUCT DATA SHEET
GENERAL SURGERY DRAPE PACK
111.001.01




ITEM 5
Reinforced Surgical Gown
Size L

ITEM 6
Hand Towels
30x40 cm



ITEM 7
OP-Tape
10x50 cm



ITEM 8
Back Table Cover
150x190 cm

(All of the product components are wrapped into back table cover)

STANDARD PRODUCT DATA SHEET
GENERAL SURGERY DRAPE PACK
111.001.01

2) DETAILS OF PRODUCT COMPONENTS

ITEM	QTY	PRODUCT	FABRIC TYPE	DIMENSION	COLOUR
1	1	Adhesive Drape <i>(60 mm x 95 cm adhesive with fingerlift)</i>	Biflex 56gr/m ² (+/-3)	150x240 cm	Blue
2	1	Adhesive Drape <i>(60 mm x 95 cm adhesive with fingerlift)</i>	Biflex 56gr/m ² (+/-3)	150x180 cm	Blue
3	2	Adhesive Drape <i>(60 mm x 95 cm adhesive with fingerlift)</i>	Biflex 56gr/m ² (+/-3)	75x90 cm	Blue
4	1	Mayo Stand Cover	Pe+nonwoven	80x145 cm	Blue
5	3	Reinforced Surgical Gown	Sms 43gr/m ² (+/-3)	Size L	Blue
6	4	Hand Towel	Cellulose	30x40 cm	White
7	1	Op-Tape	Sms+hotmelt+silicon paper	10x50 cm	Blue
8	1	Back Table Cover <i>(All of the product components are wrapped into back table cover)</i>	Biflex 56gr/m ² (+/-3)	150x190 cm	Blue

REINFORCED SURGICAL GOWN MEASUREMENT (cm)

<u>SIZE</u>	<u>A</u>	<u>B</u>	<u>C</u>	<u>D</u>	<u>E</u>	<u>F</u>	<u>G</u>
L	127	117	29	60	150	50	70

STANDARD PRODUCT DATA SHEET
GENERAL SURGERY DRAPE PACK
111.001.01

3)GENERAL INFORMATION

Product Description	GENERAL SURGERY DRAPE PACK
Reference Code	111.001.01
Packing Quantity	1 pc
Quantity Per Box	7 pcs
Sterilization	EO (Ethylene Oxide)
Certificates	EC Certificate
	ISO 13485:2016
Standards	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
Instructions For Usage	This set can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
Expiry Date	3 years

STANDARD PRODUCT DATA SHEET
GENERAL SURGERY DRAPE PACK
111.001.01

5) PACKAGING

All production	Clean Room Conditions, Medical Folding, Sterilization Bag, Peelable label,	
Labelling of product	Each unit <ul style="list-style-type: none"> - Article number, - Lot number, - Manufacturer , - Production date, - Expiry date, - EO Sterilization, - Single use, - CE, - Other customer requests, 	
Packaging	<u>Pallet system</u>	<u>Carton</u>
	Height (cm) -A	40
	Width (cm) - B	60
	Length (cm) - C	40
		
Storage Instructios	<p>HANDLING AND STORAGE: Products are sterile and must be stored in its original closed carton box, in a dry place. The products must be used as soon as the package is opened in the operation room. Do not expose to humidity and direct heat or sunlight.</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>% 10 Humidity Limitation</p> </div> </div>	