


**EXTREMITY DRAPE PACK**  
**REF.NO:152.001.01**  
**STANDARD PRODUCT DATA SHEET**  
**KY-PD-023**








STANDARD PRODUCT DATA SHEET  
EXTREMITY DRAPE PACK  
152.001.01

1) PRODUCT COMPONENTS



ITEM 1  
Extremity Drape  
200x300 cm

STANDARD PRODUCT DATA SHEET  
EXTREMITY DRAPE PACK  
152.001.01

	
<p><b>ITEM 2</b> Legging 35x80 cm (unfolded)</p>	<p><b>ITEM 3</b> Op Tape 10x50 cm</p>
	
<p><b>ITEM 4</b> Hand Towel 30x40 cm</p>	<p><b>ITEM 5</b> Mayo Stand Cover 78x145 cm (unfolded)</p>
	
<p><b>ITEM 6</b> Back Table Cover 150x190 cm <i>(All of the product components are wrapped into back table cover)</i></p>	

**STANDARD PRODUCT DATA SHEET  
EXTREMITY DRAPE PACK  
152.001.01**

**2) DETAILS OF PRODUCT COMPONENTS**

ITEM	QTY	PRODUCT	FABRIC TYPE	DIMENSION	COLOUR
1	1	Extremity Drape	Biflex 56gr/m2 (+/-3)	200x300 cm	Blue
			Extra Absorbent 110gr/m2 (+/-5)	50x50 cm	
2	1	Legging	Biflex 56gr/m2 (+/-3)	35x80 cm	Blue
3	2	Op Tape	Sms+hotmelt+silicon paper	10x50 cm	Blue
4	2	Hand Towel	Cellulose	30x40 cm	White
5	1	Mayo Stand Cover	Pe+nonwoven	78x145 cm	Blue
6	1	Back Table Cover	Biflex 56gr/m2 (+/-3)	150x190 cm	Blue


**STANDARD PRODUCT DATA SHEET  
EXTREMITY DRAPE PACK  
152.001.01**

**3)GENERAL INFORMATION**

<b>Product Description</b>	<b>EXTREMITY DRAPE PACK</b>
<b>Reference Code</b>	152.001.01
<b>Packing Quantity</b>	1 pc
<b>Quantity Per Box</b>	15 pcs
<b>Sterilization</b>	EO (Ethylene Oxide)
<b>Certificates</b>	EC Certificate
	ISO 13485:2016
<b>Standards</b>	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
<b>Instructions For Usage</b>	This set can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
<b>Expiry Date</b>	3 years

**STANDARD PRODUCT DATA SHEET**  
**EXTREMITY DRAPE PACK**  
**152.001.01**

**5) PACKAGING**

<b>All production</b>	Clean Room Conditions, Medical Folding, Sterilization Bag, Peelable label,	
<b>Labelling of product</b>	Each unit <ul style="list-style-type: none"> <li>- Article number,</li> <li>- Lot number,</li> <li>- Manufacturer ,</li> <li>- Production date,</li> <li>- Expiry date,</li> <li>- EO Sterilization,</li> <li>- Single use,</li> <li>- CE,</li> <li>- Other customer requests,</li> </ul>	
<b>Packaging</b>	<u>Pallet system</u>	<u>Carton</u>
	Height (cm) -A	40
	Width (cm) - B	60
	Length (cm) - C	40
		
<b>Storage Instructios</b>	<p><b>HANDLING AND STORAGE:</b>  Products are sterile and must be stored in its original closed carton box, in a dry place.  The products must be used as soon as the package is opened in the operation room.  Do not expose to humidity and direct heat or sunlight.</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>Humidity Limitation</p> </div> </div>	