


**ARTROSCOPY DRAPE**  
**REF.NO:250.001.01**  
**STANDARD PRODUCT DATA SHEET**  
**KY-PD-027**



**STANDARD PRODUCT DATA SHEET  
ARTHROSCOPY DRAPE  
250.001.01**

**1) PRODUCT COMPONENTS**



**ITEM 1**

Arthroscopy Drape  
200x300 cm

**2) DETAILS OF PRODUCT COMPONENTS**

ITEM	QTY	PRODUCT	FABRIC TYPE	DIMENSION	COLOUR
1	1	Arthroscopy Drape	Biflex 56gr/m2 (+/-3)	200x300 cm	Blue
		Liquid Collection Pouch	PE 45 Micron (+/-5)	75x85 cm	Transparent


**STANDARD PRODUCT DATA SHEET  
ARTHROSCOPY DRAPE  
250.001.01**

**3)GENERAL INFORMATION**

<b>Product Description</b>	<b>ARTHROSCOPY DRAPE</b>
<b>Reference Code</b>	250.001.01
<b>Packing Quantity</b>	1 pc
<b>Quantity Per Box</b>	20 pcs
<b>Sterilization</b>	EO (Ethylene Oxide)
<b>Certificates</b>	EC Certificate
	ISO 13485:2016
<b>Standards</b>	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
<b>Instructions For Usage</b>	This product can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
<b>Expiry Date</b>	3 years

**STANDARD PRODUCT DATA SHEET**  
**ARTHROSCOPY DRAPE**  
**250.001.01**

**5) PACKAGING**

<b>All production</b>	Clean Room Conditions, Medical Folding, Sterilization Bag, Peelable label,	
<b>Labelling of product</b>	Each unit <ul style="list-style-type: none"> <li>- Article number,</li> <li>- Lot number,</li> <li>- Manufacturer ,</li> <li>- Production date,</li> <li>- Expiry date,</li> <li>- EO Sterilization,</li> <li>- Single use,</li> <li>- CE,</li> <li>- Other customer requests,</li> </ul>	
<b>Packaging</b>	<u>Pallet system</u>	<u>Carton</u>
	Height (cm) -A	40
	Width (cm) - B	60
	Length (cm) - C	40
		
<b>Storage Instructios</b>	<p><b>HANDLING AND STORAGE:</b>  Products are sterile and must be stored in its original closed carton box, in a dry place.  The products must be used as soon as the package is opened in the operation room.  Do not expose to humidity and direct heat or sunlight.</p> <div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>Humidity Limitation</p> </div> </div>	

