



STANDARD PRODUCT DATASHEET-STERILE FULL PROTECTION GOWN

Reference No	81.00X.01
Product Name	Sterile Full Protection Gown
Trademark	Medcare
Manufacturer	Medcare Sağlık Ürünleri San. Ve Tic. A.Ş.


Front side of full protection gown

Back side of full protection gown



**PRODUCT DATA SHEET
FULL PROTECTION GOWN
81.00X.01**

1) PRODUCT




2) STERILE FULL PROTECTION GOWN MEASUREMENT (cm)


<u>SIZE</u>	<u>CODE</u>	<u>A</u>	<u>B</u>	<u>C</u>	<u>D</u>	<u>E</u>
S	81.001.01	110	100	28	55	140
M	81.002.01	115	105	29	58	145
L	81.003.01	127	117	29	60	150
XL	81.004.01	135	125	32	62	150
XXL	81.005.01	140	130	32	64	155
XXXL	81.006.01	145	135	34	66	160

ITEM:1

Sterile Standard Surgical Gown
Size: S-M-L-XL-XXL-XXXL



ITEM:2
Hand Towel
30x40 cm



ITEM:3
Sterilization Wrap
60x60 cm

**PRODUCT DATA SHEET
FULL PROTECTION GOWN
81.00X.01**

3) DETAILS OF PRODUCT COMPONENTS




ITEM	QTY	PRODUCT	FABRIC TYPE	DIMENSION	COLOUR
1	1	Sterile Full Protection Gown	Triplex Repellent 65gr/m ² (+/-3)	S-M-L-XL-XXL-XXXL	Blue
2	2	Hand Towel	Cellulose	30x40 cm	White
3	1	Sterilization Wrap	SS Spunbond 40gr/m ² (+/-3)	60x60 cm	Blue

4) GENERAL INFORMATION

Product Description	STERILE FULL PROTECTION GOWN
Packing Quantity	S-M-L-XL 40 pcs XXL 35 pcs XXXL 30 pcs
Sterilization	EO (Ethylene Oxide)
Certificates	EC Certificate
	ISO 13485:2016
Standards	EN ISO 14971:2019
	EN 13795-1:2019
	EN ISO 11135:2014
	EN 556-1:2001
	EN ISO 15223-1:2016
	EN ISO 11737-1:2018
	EN ISO 11737-2:2020
	EN ISO 14644-1:2015
	EN ISO 14644-2:2015
	EN 17141:2020
	EN ISO 10993-1:2020
	EN ISO 10993-5:2009
	EN ISO 10993-7:2008
	EN ISO 10993-10:2010
	EN ISO 20417:2021
	EN ISO 11607-1:2020
	EN ISO 11607-2:2020
	EN 62366-1:2015
Instructions For Usage	This product can be used during the surgery operations. This product is disposable and it prevents the bacterial transfer between the medical team and the patient, reduces the occurring of possible contamination and infection between the patient and the medical team.
Expiry Date	3 years

**PRODUCT DATA SHEET
FULL PROTECTION GOWN
81.00X.01**

5) PACKAGING

All production				
Labelling of product	Clean Room Conditions,	Medical Folding,	Sterilization Bag,	Peelable label,
Each unit				
<ul style="list-style-type: none"> - Article number, - Lot number, - Manufacturer , - Production date, - Expiry date, - EO Sterilization, - Single use, - CE, - Other customer requests, 				
Packaging				
<u>Pallet system</u>	<u>Carton</u>			
Height (cm) -A	40			
Width (cm) - B	60			
Length (cm) - C	40			
				
Storage Instructions	<p>HANDLING AND STORAGE: Products are sterile and must be stored in its original closed carton box, in a dry place. The products must be used as soon as the package is opened in the operation room. Do not expose to humidity and direct heat or sunlight.</p>			
	<div style="display: flex; justify-content: space-around; align-items: center;"> <div style="text-align: center;"> <p>45°C</p>  <p>Upper limit of temperature</p> </div> <div style="text-align: center;"> <p>% 70</p>  <p>Humidity Limitation</p> </div> </div>			